



Violet Red Bile Glucose Agar

RDM-VRBGA-01

Principle

Violet red bile glucose agar is prepared in accordance with the harmonized principles of USP/EP/IP/JP and recommended for testing bile tolerating gram-negative bacteria and microbial limit testing of pharmaceutical products and raw material used in pharmaceutical industries. Violet red bile glucose agar contains yeast extract, pancreatic digest of gelatin, bile salts, sodium chloride, glucose monohydrate, neutral red, crystal violet and agar. The pancreatic digest of gelatin serves as a source of carbon, nitrogen, vitamins and minerals. Yeast Extract supplies additional nitrogen, carbon and vitamins which stimulate bacterial growth. Glucose is a source of carbohydrate. Bile salts and crystal violet perform inhibitory action and inhibit gram positive bacteria. The fermentation of glucose to acid is indicated by the pH indicator neutral red, which changes its color to red, and by precipitation of bile acids. Organisms that rapidly attack glucose produce red purple colonies surrounded by purple haloes. Agar is a solidifying agent.

Use: Recommended for detection and enumeration of *Enterobacteriaceae*, (group of glucose fermenting) and used for testing bile tolerating gram-negative bacteria and microbial limit testing of pharmaceutical products and food products.

Contents*

Ingredients

	Gram/Litre
Yeast Extract	3.000
Pancreatic Digest of Gelatin	7.000
Bile Salts	1.500
Sodium Chloride	5.000
Glucose Monohydrate	10.000
Neutral Red	0.030
Crystal Violet	0.002
Agar	15.000
pH at 25°C	7.4 ±0.2

* Formula adjusted for optimum performance and parameters

Directions: Dissolve 41.5 grams in 1000 ml distilled water. Heat to boil for 2 minutes to dissolve ingredients completely, **Do Not Autoclave**. In aseptic condition immediately cool it to 42-45°C and distribute in sterile petri plates and allow solidify. Ensure complete solidification and inoculate test sample.

Specimens types analyzed

Pharmaceutical samples, clinical and non-clinical samples etc.

Precautions to be taken

These microbial media are intended for the in-vitro use only. All the handling, experiments, storage, and discarding should be performed with the help of skilled and knowledgeable technicians and as per the established guidelines. The material should be disposed only after proper sterilization by autoclaving. Please go through the MSDS of the media to avoid any accidents or in emergency.

Performance and Evaluation

The expected performance of the medium is liable to use as per the direction on the label when stored at optimum conditions and within expiry date.

Quality Control

Appearance	Pink-beige colored free flowing, homogeneous powder
Reaction of 4.15% solution	7.4 ±0.2 at 25 °C
pH	7.20- 7.60
Gelling	Firm comparable with 1.5% agar gel
Color and clarity of ready medium	Reddish-purple, slightly opalescent gel
Growth Promotion properties	Best at ≤ 100 CFU at 32-37 °C for 18-72 h
Indicative properties	Optimum at ≤ 100 CFU at 32-37 °C for 18-48 h
Negative control	Performed using sterile distilled water

Different Microbial Response

Organism	Inoculum CFU	Growth	Recovery	Colony color	Incubation period
<i>Escherichia coli</i> (ATCC 8739)	50-100	Luxurious	60-70%	Pink-red with bile precipitate	30-37°C, 18-24 hours
<i>Pseudomonas aeruginosa</i> (ATCC 27853)	50-100	Luxurious	60-70%	Pink to purple	30-37°C, 18-24 hours
<i>Staphylococcus aureus</i> (ATCC 25923)	50-100	Inhibited	---	--	30-37°C, 18-24 hours

Storage and Shelf Life

Hygroscopic; keep container tightly closed. Store in cool dry place.

Disposal: To avoid the contamination or propagation of any hazardous microbes the used, unusable or modified preparation of this product must be disposed after autoclaving after completion of task.

Reference

1. Atlas, R. M. (2005). *Handbook of media for environmental microbiology*. CRC press.
2. *British Pharmacopoeia*, (2011), The Stationery office British Pharmacopoeia
3. *Difco Manual* (1998). 11th Edition. Difco Laboratories., Division of Becton Dickinson and Company, Sparks, Maryland, USA.
4. *European Pharmacopoeia*, (2011), European Dept. for the quality of Medicines.
5. *Indian Pharmacopoeia*, (2018), Govt. of India, the Controller of Publication, New Delh
6. Rand, M. C., Arnold E. Greenberg, and Michael J. Taras, (1976), *Standard methods for the examination of water and wastewater*. Prepared and published jointly by American Public Health Association, American Water Works Association, and Water Pollution Control Federation.
7. *The Japanese Pharmacopoeia*, 17th Ed. (2016), The Ministry of Health, Labour And Welfare
8. *The United States Pharmacopoeia*, (2014), The United States Pharmacopeial Convention. 12601 Twinbrook Parkway, Rockvukke, MD 20852.

Disclaimer

The information contained in the technical data sheet is to the best of our knowledge is accurate and true based on the research and development work carried out by **ReadyMED**[®], Chaitanya Agro Biotech, Malkapur, Maharashtra. The products are neither intended for any therapeutic use for animal or human nor for any other *in-vivo* applications. The **ReadyMED**[®] products are only meant to be used for the laboratory, diagnostic, research, or further manufacturing purpose only. These technical outcomes should not be considered as the warranty of any kind expressed or implied, and no liability is accepted for infringement of any patent.

CHAITANYA AGRO BIOTECH PVT. LTD. An ISO 11134:2014, ISO 13485:2016, ISO 9001:2015 CE, CIN NO.: U24210MH1995PTC095220S, S. No. 120/2, Laxmi Nagar, Umbarnala Road, Malkapur-443101, Dist.: Buldana (M.S.) India. Customer Care +91- 8669083859
rdmsales@chaitanyagroupindia.com, mkt.cabt@chaitanyagroupindia.com, www.chaitanyagroupindia.com