



Eosine Methylene Blue Broth

RDM-EMBB-01

Principle

Eosine Methylene Blue Broth composed of peptone, lactose, sucrose, dipotassium phosphate, Eosin and methylene blue. Peptone as a source of carbon, nitrogen, vitamins and minerals. Lactose and sucrose are fermentable carbohydrates. Dipotassium phosphate acts as buffer. Eosin and methylene blue are dyes provide selectivity to media and act as pH indicator and inhibitors. The media also used for cultivation of lactose fermenters and non-lactose fermenters. Some gram-negative bacteria can not able to metabolize lactose or metabolize lactose very slowly. Hence, the sucrose is used as secondary carbohydrate source. While the lactose non-fermenters such as *Salmonella sp.* raise the pH of surrounding by deamination of protein thus resulting in change of pH, which solubilizes the methylene blue-eosin complex resulting, color of media changes from light purple to pink or colorless. While the coliforms or lactose fermenters produce acid and decrease pH and media color changes form purple to dark purple with greenish sheen, such as *Escherichia coli* forms characteristic green colored metallic sheen. Gram-positive bacteria are partially inhibited on the medium.

Use: Recommended for the differential and cultivation of the gram-negative Enteric Bacilli from clinical and non-clinical samples.

Contents*

Ingredients

	Gram/Litre
Peptone	10.000
Lactose	5.000
Sucrose	5.000
Dipotassium Phosphate	2.000
Eosin Y	0.400
Methylene Blue	0.065
pH at 25°C	7.2 ±0.2

* Formula adjusted for optimum performance and parameters

Directions: Dissolve 22.50 grams in 1000 ml distilled water check. Boil to dissolve the medium completely and distribute aseptically. Sterilize by autoclaving at 15 lbs pressure (121 °C) for 15 min, cool it to 42-45 °C and inoculate test sample aseptically.

Specimens types analyzed

The water samples, soil samples, food and dairy samples, pharmaceutical samples, clinical samples etc.

Precautions to be taken

These microbial media are intended for the in-vitro use only. All the handling, experiments, storage, and discarding should be performed with the help of skilled and knowledgeable technicians and as per the established guidelines. The material should be disposed only after proper sterilization by autoclaving. Please go through the MSDS of the media to avoid any accidents or in emergency.

Performance and Evaluation

The expected performance of the medium is liable to use as per the direction on the label when stored at optimum conditions and within expiry date.

Quality Control

Appearance	Light purple colored free flowing, homogeneous powder
Reaction of 2.25% solution	7.2 ±0.2 at 25 °C
pH	7.00- 7.40
Color and clarity of ready medium	Light purple colored opalescent solution
Growth Promotion properties	Best at ≤ 100 CFU at 32-37 °C for 18-72 h
Indicative properties	Optimum at ≤ 100 CFU at 32-37 °C for 18-48 h
Negative control	Performed using sterile distilled water

Different Microbial Response

Organism	ATCC	Inoculum	Growth	Colony Color	Incubation
<i>Escherichia coli</i>	8739	50-100	Luxurious	Dark purple with green ting	33-37°C , 18-48 h
<i>Salmonella typhi</i>	14028	50-100	Good	Colorless	33-37°C , 18-48 h

Storage and Shelf Life

Hygroscopic; keep container tightly closed. Store in cool dry place.

Disposal: To avoid the contamination or propagation of any hazardous microbes the used, unusable or modified preparation of this product must be disposed after autoclaving after completion of task.

Reference

1. Atlas, R. M. (2004). *Handbook of microbiological media*. CRC press.
2. Atlas, R. M. (2005). *Handbook of media for environmental microbiology*. CRC press.
3. *Difco Manual* (1998). 11th Edition. Difco Laboratories., Division of Becton Dickinson and Company, Sparks, Maryland, USA.
4. Rand, M. C., Arnold E. Greenberg, and Michael J. Taras. *Standard methods for the examination of water and wastewater*. Prepared and published jointly by American Public Health Association, American Water Works Association, and Water Pollution Control Federation., 1976.

Disclaimer

The information contained in the technical data sheet is to the best of our knowledge is accurate and true based on the research and development work carried out by **ReadyMED**[®], Chaitanya Agro Biotech, Malkapur, Maharashtra. The products are neither intended for any therapeutic use for animal or human nor for any other *in-vivo* applications. The **ReadyMED**[®] products are only meant to be used for the laboratory, diagnostic, research, or further manufacturing purpose only. These technical outcomes should not be considered as the warranty of any kind expressed or implied, and no liability is accepted for infringement of any patent.

CHAITANYA AGRO BIOTECH PVT. LTD. An ISO 11134:2014, ISO 13485:2016, ISO 9001:2015 CE , CIN NO.: U24210MH1995PTC095220S,
S. No. 120/2, Laxmi Nagar, Umbarnala Road, Malkapur-443101, Dist.: Buldana (M.S.) India. Customer Care +91-8669083859
rdmsales@chaitanyagroupindia.com, mkt.cabt@chaitanyagroupindia.com, www.chaitanyagroupindia.com