



Lowenstein Jensen Medium Base

RDM-LJMB-01

Principle

Lowenstein-Jensen medium base, is used for isolation of acid-fast microorganisms like Mycobacterium species. The base medium is composed of L-asparagine, monopotassium phosphate, magnesium sulphate, magnesium citrate, potato starch soluble and malachite green. L-asparagine is source of nitrogen. Monopotassium phosphate and magnesium sulfate enhance organism growth and act as buffers. Potato starch is source of carbon and absorb toxic metabolites, malachite green is used to prevent the growth of most contaminants and encourages the growth of mycobacteria.

For better result the media can be fortified with egg suspension, glycerol, sodium chloride, penicillin and nalidixic acid alone or in combination, Egg suspension, which provides fatty acids and protein required for the metabolism of mycobacteria. Glycerol serves as carbon source. Sodium chloride helps to differentiate rapid-growing mycobacteria from slow growers. Penicillin and nalidixic acid are added to decrease contamination.

Use: Recommended for isolation, cultivation and differentiation of mycobacteria.

Contents*

Ingredients

	Gram/Litre
L-Asparagine	3.60
Monopotassium Phosphate	2.40
Magnesium Sulphate	0.24
Magnesium Citrate	0.60
Potato Starch Soluble	30.00
Malachite Green	0.40

* Formula adjusted for optimum performance and parameters

Directions: Dissolve 37.3 gm of the medium in 600 ml of distilled water containing 12 ml of glycerol. Heat the medium to dissolve completely, autoclave at 121°C for 15 minutes. Prepare 1000 ml of a uniform suspension of fresh eggs under aseptic conditions. Swirl gently to avoid introducing air into the suspension. Aseptically mix the 1000 ml of egg suspension with 600 ml of the sterile Lowenstein-Jensen Medium base cooled to 50 – 60°C, avoiding air bubbles. Dispense the finished medium into sterile screw-cap test tubes. Place the tubes in a slanted position and heat at 85°C for 45 minutes to coagulate the medium.

Specimens types analyzed

Pharmaceutical samples, clinical and non-clinical samples etc.

Precautions to be taken

These microbial media are intended for the in-vitro use only. All the handling, experiments, storage, and discarding should be performed with the help of skilled and knowledgeable technicians and as per the established guidelines. The material should be disposed only after proper sterilization by autoclaving. Please go through the MSDS of the media to avoid any accidents or in emergency.

Performance and Evaluation

The expected performance of the medium is liable to use as per the direction on the label when stored at optimum conditions and within expiry date.

Quality Control

Appearance	A fine, free flowing, hygroscopic dried powder. blue to greenish blue powder
Color and clarity of ready medium (after addition of egg suspension)	The mixture of sterile basal medium and whole egg emulsion, coagulates to yield pale bluish green colored, opaque smooth slant
Growth Promotion properties	Best at ≤ 100 CFU at 32-37 °C for 18-72 h
Indicative properties	Optimum at ≤ 100 CFU at 32-37 °C for 18-48 h
Negative control	Performed using sterile distilled water

Different Microbial Response

Organism	Inoculum CFU	Growth	Colony characteristic	Incubation
<i>Mycobacterium avium</i> (ATCC 25291)	50-100	Luxuriant	smooth, non-pigmented colonies	33-37°C 2-4 weeks in presence of CO ₂

Storage and Shelf Life

Hygroscopic; keep container tightly closed. Store in cool dry place.

Disposal: To avoid the contamination or propagation of any hazardous microbes the used, unusable or modified preparation of this product must be disposed after autoclaving after completion of task.

Reference

1. Atlas, R. M. (2005). *Handbook of media for environmental microbiology*. CRC press.
2. *Difco Manual* (1998). 11th Edition. Difco Laboratories., Division of Becton Dickinson and Company, Sparks, Maryland, USA.

Disclaimer

The information contained in the technical data sheet is to the best of our knowledge is accurate and true based on the research and development work carried out by **ReadyMED®**, Chaitanya Agro Biotech, Malkapur, Maharashtra. The products are neither intended for any therapeutic use for animal or human nor for any other *in-vivo* applications. The **ReadyMED®** products are only meant to be used for the laboratory, diagnostic, research, or further manufacturing purpose only. These technical outcomes should not be considered as the warranty of any kind expressed or implied, and no liability is accepted for infringement of any patent.

CHAITANYA AGRO BIOTECH PVT. LTD. An ISO 11134:2014, ISO 13485:2016, ISO 9001:2015 CE, CIN NO.: U24210MH1995PTC095220S,
S. No. 120/2, Laxmi Nagar, Umbarnala Road, Malkapur-443101, Dist.: Buldana (M.S.) India. Customer Care +91-8669083859
rdmsales@chaitanyagroupindia.com, mkt.cabt@chaitanyagroupindia.com, www.chaitanyagroupindia.com