



Technical Data Sheet

Format no: RDM/TD/P-03

Revision: 01 / 2024

Cas Amino Acid Technical Grade

RDM-P-03

Principle: Cas amino acids, is complete hydrolyzed product of casein containing minute concentration of sodium chloride and iron. It is less refined product of cas amino acid purified. In media preparations is used as replacement of pure amino acids and small peptides. It has been utilized in diverse media compositions used for cultivation of lactobacillus and parasite (*Entamoeba histolytica*) and nematodes medium.

Use: Recommended to use as culture media ingredient in variety of media preparations used for production of enzymes, steroids, vaccines and other pharmaceutical, agricultural economical products.

Quality Control

Appearance	Off white colored homogeneous free flowing, powder
Solubility (2%)	Soluble in distilled water
Clarity	Pale yellow colored clear solution without haziness at 2 % concentration
pH	6.00 – 7.50 at 25°C
Loss on drying	NMT 7.00% as estimated by AOAC method

Chemical analysis

Total Nitrogen	NLT 8.00%
Amino Nitrogen	NLT 4.00%
Residue on ignition	NMT 39.00%

Bacteriological testing: Bacteriological tests are carried out as per USP 32, NF26 where respective medium is prepared by using Bi-peptone under test.

Test for pathogens:

Total Plate Count	NMT 10000 cfu per gram.
Yeast & Molds	Absent per 10 gram.
<i>Escherichia coli</i>	Absent per 10 gram.
<i>Salmonella</i>	Absent per 10 gram.
<i>Staphylococcus aureus</i>	Absent per 10 gram.

Culture response: Cultural response observed after incubation at 35-37°C for 24 hours by using 2% Cas amino acid, 0.5% sodium chloride and 1.5% agar, vitamin mixture 0.001% in water, pH 7.2-7.4.

<i>Escherichia coli</i> (ATCC 25922)	Luxurious growth
<i>Salmonella typhimurium</i> (ATCC 14028)	Luxurious growth
<i>Pseudomonas aeruginosa</i> (ATCC 10145)	Luxurious growth

Storage and Shelf Life

Store below 30°C in tightly sealed jar or container. Use before expiry date on the label. Expected performance when stored at optimum conditions and within expiry date.

Disposal: To avoid the contamination or propagation of any hazardous microbes used, unusable or modified preparation of this product must be disposed after autoclaving or incineration after completion of task.

Disclaimer

The information contained in the technical data sheet is to the best of our knowledge is accurate and true based on the research and development work carried out by **ReadyMED®**, Chaitanya Agro Biotech, Malkapur, Maharashtra. The products are neither intended for any therapeutic use for animal or human nor for any other *in-vivo* applications. The **ReadyMED®** products are only meant to be used for the laboratory, diagnostic, research, or further manufacturing purpose only. These technical outcomes should not be considered as the warranty of any kind expressed or implied, and no liability is accepted for infringement of any patent.

CHAITANYA AGRO BIOTECH PVT. LTD.

An ISO 11133:2014, ISO 13485:2016, ISO 9001:2015 ॐ , CIN NO.: U24210MH1995PTC095220S,

S. No. 120/2, Laxmi Nagar, Umbarnala Road, Malkapur-443101, Dist.: Buldana (M.S.) India. Customer Care +91- 8669083859

rdmsales@chaitanyagroupindia.com, mkt.cabt@chaitanyagroupindia.com, www.chaitanyagroupindia.com