



Technical Data Sheet

Format no: RDM/TD/H-TSA-01

Revision: 01 / 2024

Tryptic Soya Agar (Harmonized)

RDM-H-TSA-01

Principle

Tryptic Soya Agar (Harmonized) is used for cultivation of a wide variety of organisms. The medium is a highly nutritious composed of tryptone, soya peptone, sodium chloride, dextrose, dipotassium hydrogen phosphate and agar. Tryptone, papain digest of soy peptone provides ample amount of the carbon and nitrogen source and trace. Dextrose is carbon source. Phosphate is buffering agent. Sodium chloride maintains the osmotic balance and agar is used as solidifying agent. Media is specially used for the cultivation of *Salmonella typhimurium* from clinical and non-clinical samples, food and dairy products. In addition to that it is also used for isolation, cultivation and maintenance of other fastidious organisms too.

Use: Recommended for the isolation and cultivation and maintenance of *Salmonella typhimurium* from clinical and non-clinical samples.

Contents*

Ingredients	Gram/Liter
Tryptone	17.00
Soya Peptone	3.00
Sodium Chloride	5.00
Dextrose	2.50
Dipotassium hydrogen phosphate	2.50
Agar	15.00
pH at 25°C	7.3 ±0.2

* Formula adjusted for optimum performance and parameters

Directions: Dissolve 45.00 grams in 1000 ml distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121 °C) for 15 min, cool it to 42-45 °C and distribute aseptically in petri plates and allow to solidify. Ensure complete solidification and inoculate test sample aseptically.

Specimens types analyzed

Pharmaceutical samples, clinical and non-clinical sample.

Precautions to be taken

All the handling, experiments, storage, and discarding should be performed with the help of skilled and knowledgeable technicians and as per the established guidelines. The material should be disposed only after proper sterilization by autoclaving. Please go through the MSDS of the media to avoid any accidents or in emergency.

Performance and Evaluation

The expected performance of the medium is liable to use as per the direction on the label when stored at optimum conditions and within expiry date.

Quality Control

Appearance	Light beige colored, free-flowing, homogeneous
Reaction of 4.5% solution	7.30 ±0.2 at 25 °C
pH	7.10- 7.50
Gelling	Firm comparable with 1.5% agar gel
Color and clarity of ready medium	Light amber, clear opalescent gel
Growth Promotion properties	Best at ≤ 100 CFU at 32-37 °C for 18-72 h
Indicative properties	Optimum at ≤ 100 CFU at 32-37 °C for 18-48 h
Negative control	Performed using sterile distilled water

Different Microbial Response: Prepare media as per the label directions. Inoculate and incubate the plates at 35 ± 2.0°C for 24-48 hours.

Organism	ATCC	Growth	Recovery
<i>Salmonella typhimurium</i>	14028	Luxuriant	≥ 70%
<i>Salmonella abony</i>	NCTC-6017	Luxuriant	≥ 70%
<i>Salmonella enterica</i>	9150	Luxuriant	≥ 70%
<i>Salmonella enteritidis</i>	13076	Luxuriant	≥ 70%

Storage and Shelf Life: The product is highly hygroscopic; keep the container tightly closed at all times and store it properly as per the conditions mentioned on the label. The declared expiry is valid only when stored as per the conditions mentioned on the label.

Note: Sterilize media immediately after reconstitution.

Disposal: To avoid the contamination or propagation of any hazardous microbes the used, unusable or modified preparation of this product must be disposed after autoclaving after completion of task.

Reference

1. Atlas, R. M. (2005). *Handbook of media for environmental microbiology*. CRC press.
2. *Difco Manual* (1998). 11th Edition. Difco Laboratories., Division of Becton Dickinson and Company, Sparks, Maryland, USA.

Disclaimer

The information contained in the technical data sheet is to the best of our knowledge is accurate and true based on the research and development work carried out by **ReadyMED®**, Chaitanya Agro Biotech, Malkapur, Maharashtra. The products are neither intended for any therapeutic use for animal or human nor for any other *in-vivo* applications. The **ReadyMED®** products are only meant to be used for the laboratory, diagnostic, research, or further manufacturing purpose only. These technical outcomes should not be considered as the warranty of any kind expressed or implied, and no liability is accepted for infringement of any patent.

CHAITANYA AGRO BIOTECH PVT. LTD.

An ISO 11133:2014, ISO 13485:2016, ISO 9001:2015  , CIN NO.: U24210MH1995PTC095220S,
S. No. 120/2, Laxmi Nagar, Umbarnala Road, Malkapur-443101, Dist.: Buldana (M.S.) India. Customer Care +91-8669083859
rdmsales@chaitanyagroupindia.com, mkt.cabt@chaitanyagroupindia.com, www.chaitanyagroupindia.com